


Automated, non-invasive culture, and evaluation system for iPS cells under neural differentiation process

Chieko Nakada¹, Momotaro Ishikawa¹, Hiroaki Kii¹, Takayuki Uozumi¹, Hirotsada Watanabe¹, Yoichi Yamazaki², Takayuki Kitogo³, Yasuyuki Asai³, Yasujiro Kiyota¹
¹ Instruments company, NIKON CORPORATION, Yokohama, Japan, ² NIKON SYSTEMS, Inc., Tokyo, Japan, ³ ReproCELL, Inc., Kanagawa, Japan


Introduction

Differentiation process of human iPS cells toward DA neuron



- It takes more than ~4 weeks in general to induce neurons from iPS/ES cells.
- "Neural rosette" and "neurites" are the two morphological checkpoints in the neural differentiation process.
- We aimed to develop an integrated culture, observation and evaluation system which enables quantitative and label-free analysis of cells being differentiated from iPS cells into neurons under incubation condition.
- New image processing methods were developed to detect neural rosette structures and neurons (cell bodies and neurites);
- The developed methods were then tested and applied to estimate the differentiation status.

Methods



Schematic drawing showing experimental workflow.

Cell culture

- iPS cells were differentiated from different lines of iPS cells into Dopaminergic neurons.
- Cells at rosette stage were plated on 6 well-plates. To culture neural precursors for further differentiation and maturation, rosette-stage cells were collected, dissociated, and plated on 24-well plate at densities of 15 and 7.5 x 10⁴ cells/well. Cells were kept in BioStation CT (NIKON CORPORATION) at 37°C, 5% CO₂.

Time-lapse image acquisition

- All imaging was performed using BioStation CT (NIKON CORPORATION) at 37°C, 5% CO₂.
- For rosette-stage cells, phase contrast images were captured at a magnification of 2x. To obtain high contrast at any height of the cell structure, z-stack of 5μm x 40 images was acquired at each x-y position. These z-stack images were used in the image analysis to create EDF (extended depth of focus) image. 25 x-y positions were observed at each well, so that the squared area of 19.2 x 19.2 mm was covered. Images were captured every 6 hours for 5 to 7 days.

- For observation of neurons and precursors, phase contrast images were obtained at 10x and 20x. Neurite outgrowth was observed by capturing images at every 30 min or 60 min for 5 days for 10x and 20x observation, respectively. After time lapse observations, neurons were fixed, stained with a monoclonal antibody against MAP2 (a neurite marker protein) followed by Alexa 488-labeled secondary antibodies, and the fluorescence signals were observed with the GFP channel.

Image analysis

- All image analysis was performed using functions built-in software (CL-Quant, NIKON CORPORATION) [1].




BioStation CT


Results 1.

Neural rosette structure detection and analysis

Algorithm Development




Algorithm Validation



Validation of the developed algorithm. Left, correlation in the number of rosette-like structures between human eye (manual count, n=3) and the algorithm. It showed a high correlation coefficient of 0.98. Right, Bland-Altman plot.


Application: Rosette-forming efficiency: lot-to-lot variation




Results 2.

Neurite Detection and Outgrowth Analysis


Algorithm Development



Detection Results




Detection Results



Neurite elongation was monitored over time (30 min x 5 days). Two samples (derived from different rosette samples (lot A and B) were compared and there was no significant difference throughout the observation period. The top bar graphs show the values at the final time point.

Application: Monitoring cell survival over time



Reference

1. Alwrth, SV, Watanabe, H. & Lee, JSJ. Teachable, High-Content Analytics for Live-Cell, Phase Contrast Movies. J Biomol Screening 8, 968-977, 2010.

Conclusion

- Phase contrast, live cell, Image-based analysis was performed with new algorithms developed for neural rosette-structure and neurite detection.
- These algorithms detected rosettes in high correlation with human eye.
- The neurites detected based on the phase-contrast images highly matched with those based on fluorescent signals from neurite marker.
- These analysis would become useful tools for variety of applications, such as estimation of differentiation efficiency in different protocols, toxicity testing, and benchmarks for quality controls.